

## FLOCARE® ACCESSORIES

The Flocare® accessories range includes feeding connectors, adaptors, containers, and pump accessories to aid enteral feeding.

All Flocare® Feeding Tubes, Flocare® Feeding Sets and ENFit™ enteral syringes are ENFit™ compatible.

ADAPTORS	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
Flocare® Transition Step Connector	Designed to be attached to an ENFit <sup>™</sup> feeding set to enable connection to a non-ENFit <sup>™</sup> feeding tube. Recommended for 24-hour, single use, single patient only.	89828	30
Flocare® ENFit™ Bolus Adaptor	Designed to be inserted into ENPlus enteral feeding packs to enable connection to an ENFit™ syringe to allow for bolus feeding. Recommended for 24-hour, single tube feed pack use.	89740	30
Flocare® Universal Pack Adaptor	Designed to be used alongside Flocare® Infinity Pack feedings sets to enable connection to an enteral bottle or Flocare® container. Recommended for 24-hour, single use, single patient only.	70063	30

CONTAINERS	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
	Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Recommended for 24 hour, single use only.	40455	10
Flocare® Container 1L			

MOBILE ACCESSORIES	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
Flocare® Infinity GO Frame	Table-top, lightweight, plastic frame that easily fixes to the Flocare® Infinity pump and will allow enteral pack / bottle, Op Tri bottles or Flocare® containers to be attached, for ease of delivery.	43247	1
	Designed to optimise patient mobility and product protection when tube feeding. Capacity to carry both the Flocare® Infinity pump, enteral feed and GO frame at the same time.	44497	1
Flocare® Infinity GO Bag	at the same time.		

For more information, please call: Nutricia Clinical Care Line 0800 636 228

